



Focus on
Incubators

LINAK.COM/MEDLINE-CARELINE


LINAK® 
WE IMPROVE YOUR LIFE

Flexible functionality and ergonomics for incubators

Care and nurturing of premature born babies call for intelligent and reliable solutions. **Electric adjustable incubators** provide improved comfort and ergonomics for work staff and relatives, as the height can be adjusted for individual use. With an electric actuator solution using the OpenBus™ technology from LINAK, there are several ways to **increase the value of Incubators.**

LINAK offers:

- **High-quality electric solutions**
- **Add-on accessories via OpenBus™ technology**
- **Flexible solutions according to your needs**




Facts

LINAK has more than 30 years of experience within the actuator business and can provide electric actuator solutions of high quality that enable you to create a unique and differentiated product via our OpenBus technology.

Our solutions have been developed and tested as systems. This ensures the functioning and the safety of our systems.

LINAK offers intelligent system solutions via OpenBus™



The diagram shows a central incubator unit with a blue control panel. To the left is an OpenBus control box. Below the incubator is a foot switch. To the right are three BL1 columns of different heights. Above the incubator is an SCO scale. A control panel with a display and buttons is also shown to the right of the incubator.

OpenBus control boxes

- Add-on accessories possible
- Low weight
- Cable management

FS3 Foot switch

- Intuitive - easy to use
- Aesthetic design and expandable modules
- Cleaning-friendly

SCO scale

- Fast and easy weighing in the incubator
- Add-on accessory
- Customised front covers possible

BL1 columns

- Aluminium column
- Low noise level
- Up to 2000 N in push

Optional cables inside the column for aesthetic design, high safety and cleaning-friendly

OpenBus™ accessories





GLOBAL PRESENCE

LINAK® has a well-developed sales and service organisation in Europe, Americas, Asia and Australia. Therefore, we can assist you and your customers locally, under the global sales concept idea:

Be global, act local



For further information, please visit our website:

www.linak.com/medline-careline

Please scan the QR code for direct access to Incubators

TERMS OF USE

The user is responsible for determining the suitability of LINAK products for a specific application. LINAK takes great care in providing accurate and up-to-date information on its products. However, due to continuous development in order to improve its products, LINAK products are subject to frequent modifications and changes without prior notice. Therefore, LINAK cannot guarantee the correct and actual status of said information on its products. While LINAK uses its best efforts to fulfil orders, LINAK cannot, for the same reasons as mentioned above, guarantee the availability of any particular product. Therefore, LINAK reserves the right to discontinue the sale of any product displayed on its website or listed in its catalogues or other written material drawn up by LINAK.

All sales are subject to the Standard Terms of Sale and Delivery for LINAK. For a copy hereof, please contact LINAK.

LINAK has a world-class sales and service organisation.
Today we are present in 35 countries all over the world.
For further information, please visit our website: www.linak.com

LINAK® 
WE IMPROVE YOUR LIFE